










The Design								
The Product Line	<b>The Red Line</b>							
One-Knob Data Entry								
Diameter of the Chamber in mm	400	400	400	400	520	520	400	520
Model	HiCen T	HiCen TR	HiCen F	HiCen FR	HiCen GT	HiCen GR	HiCen SR	HiCen XL
Table / Floor	Table top	Table top	low Floor Model	low Floor Model	Table top	Floor Model	Floor Model	Floor Model
Max. Speed	14.000	16.000	14.000	18.000	14.000	18.000	21.000	21.000
g-Force	21.913	26.790	21.913	34.779	21.913	33.688	50.743	50.743
Refrigerated	no	yes	no	yes	no	yes	yes	yes
Max. Volume Swing-out in ml	4 x 500	4 x 500	4 x 500	4 x 500	4 x 1000	4 x 1000	4 x 500	4 x 1000
Max. Volume Fixed-Angle in ml	6 x 250	6 x 250	6 x 250	6 x 250	4 x 1000	4 x 1000	4 x 1.000	6 x 1.000
Max. Capacity 1,5/2,2 ml Tubes	72	72	72	72	72	72	72	72
Max. Capacity 15 ml Falcon® Tubes	56	56	56	32	76	76	56	76
Max. Capacity 50 ml Falcon® Tubes	28	28	28	28	32	32	28	32

Herolab reserves the right to change technical specification without prior notice. 04/2016

