

9. Rotor Matching Guide for Tubes and Bottles

Herolab

Tubes and Bottles for Highspeed Rotors

Material Abbreviations
 B - Borex Glass® PA - Polyallomer PF - Polyfo®
 C - Corex Glass® PC - Polycarbonate PP - Polypropylene
 D - Duran Glass® PE - Polyethylene PS - Polystyrene

NOMINAL CAPACITY ROTOR MATCHING GUIDE *For Commonly Used Rotors*

10 ml Herolab AFH 20.180 Beckman JA-21 Sigma 12157 12111	14 ml Sovall ST-12 Kontron A1224	16 ml Beckman JA-201 Sovall SW-24 Kontron A 24 24	30 ml Sigma 12158	50 ml Herolab AFH 8.180 Kontron A 8.24 Sigma 12156 12150 Beckman JA-20 JS-13 JA-17 JS-75	80 ml Kontron A 8.20 AS 4.13-2 AV 8.19 Beckman JA-18 JS-13	250 ml Herolab AFH 4.080 Sigma 12159 12155	250 ml Herolab AFH 6.180 Kontron AS 4.13-2 AS 4.7 Sigma 12159 12155	500 ml Kontron A 6.9 Beckman JA-10 Sovall GS-3
---	---	---	-----------------------------	---	--	--	---	--

A = Adaptor required: See individual rotor data books or the Herolab Centrifuge Labware Catalog.

10 ml 16 x 80 mm Cat. No. PA - 25 24 11 24,000 PC - 25 24 01 24,000 PE - 25 24 51 20,000 PP - 25 24 21 24,000	13 ml 18 x 102 mm Cat. No. B - 25 11 00 A	14 ml 19 x 75 mm Cat. No. PA - 25 26 90 24,000 PC - 25 25 70 24,000 PE - 25 25 80 20,000 PP - 25 25 80 20,000	15 ml 18 x 102 mm Cat. No. C - 25 10 90 9,000/12,500 A	16 ml 18 x 100 mm Cat. No. PA - 25 27 20 24,000 PC - 25 27 00 24,000 PE - 25 27 50 20,000 PP - 25 27 10 20,000	30 ml 24 x 105 mm Cat. No. B - 25 11 10 C - 25 11 40 A	45 ml 28 x 102 mm Cat. No. B - 25 12 60 A	50 ml 29 x 103 mm Cat. No. PA - 25 31 60 24,000 PC - 25 31 30 20,000 PE - 25 31 51 20,000 PP - 25 31 10 20,000	50 ml 29 x 103 mm Cat. No. PA - 25 29 11 24,000 PC - 25 29 01 12,000 PE - 25 29 21 15,000 PP - 25 29 05 24,000	80 ml 38 x 105 mm Cat. No. PA - 25 32 45 20,000 PC - 25 32 41 12,000 PE - 25 32 47 15,000 PP - 25 32 49 15,000 PP - 25 32 43 20,000
--	--	--	---	---	--	--	---	---	---

250 ml 62 x 122 mm Cat. No. PA - 25 34 21 14,400 PC - 25 34 01 14,400 PE - 25 34 31 13,500 PF - 25 34 39 14,400 PP - 25 34 11 14,400 PP - 25 34 12/Black 14,400	250 ml 62 x 120 mm Cat. No. PA - 25 34 45 14,400 PC - 25 34 41 14,400 PE - 25 34 47 13,500 PF - 25 34 49 12,000 PP - 25 34 43 14,400	290 ml 62 x 137 mm Cat. No. PA - 25 34 85 14,400 PC - 25 34 81 14,400 PE - 25 34 87 13,500 PP - 25 34 83 14,400	500 ml 70 x 165 mm Support Ring Cat. No. 25 52 41 Cat. No. PA - 25 35 21 9,000 PC - 25 35 01 9,000 PE - 25 35 31 9,000 PP - 25 35 12 / Black 9,000
--	--	--	--



Herolab GmbH
 Ludwig-Wagner-Str.12
 D-69168 Wiesloch
 Tel-Nr.: +49 (0) 6222/5802-0
 Fax-Nr.: +49 (0) 6222/5802-34
 E-Mail: herolab_gmbh@t-online.de
 Internet: http://www.Herolab.de

The information in this wall chart is intended as a general guide only. Since the performance of centrifuge labware depends on many different variables, no absolute guarantee of performance can be given. Should there be any doubts, users should experiment with the centrifuge labware under their own operating conditions prior to the actual test. Specifications are subject to change.

All products built by Herolab carry the CE sign. Herolab has its own Quality Management System which is conform with ISO 9001. All prices without VAT, ex works. – Herolab reserves the right to change technical specifications and prices without prior notice.