

Herolab

Tubes and Bottles for Ultraspeed Rotors

Made in Switzerland
Quality assured



2ml
8x51mm

PA, thin ③ 25 20 70
PC, thick 25 20 20



4,4ml
11x60mm

PA, thin ③ 25 21 40
PC, thick 25 21 50
PC, thick, thr 25 22 20/A①
Polyclear 25 21 80



5ml
13x51mm

PA, thin 25 21 90
PC, thick 25 22 60
PC, thick 25 22 40
Polyclear 25 22 50



6,5ml
13x64mm

PA, thin ②③ 25 23 40
PA, thick 25 23 50
PC, thick 25 22 80
PC, thick, thr 25 22 90 / A①
Polyclear 25 23 80



13,2ml
14x89mm

PA, thin 25 23 90



13,5ml
16x76mm

PA, thin ②③ 25 24 30
PC, thick ② 25 24 70
PC, thick ② 25 24 80
PC, thick, thr 25 24 90 / A①
Polyclear 25 24 50



14ml
14x95mm

PA, thin 25 23 60
PC, thick 25 25 60
Polyclear 25 25 30




17ml
16x102mm

PA, thin ③ 25 26 20
PC, thick ② 25 26 40
Polyclear 25 27 60




38,5ml
25x89mm

PA, thin ②③ 25 30 50
PA, thick ② 25 30 60
PC, thick ② 25 30 40
PC, thick, thr 25 30 20 / A①
Polyclear 25 30 70



94ml
38x102mm

PC, thick ② 25 32 80
PC, thick, thr 25 32 90 / A①
Polyclear 25 32 70



250ml
60x122mm

PA 25 34 25

Material Abbreviations

PA - Polyallomer
PC - Polycarbonate
thin - thinwall
thick - thickwall
thr - threaded
A - Assembly / Tube & Screw Cap

Caution:
Obey rotor manufacturers
speed limits!

Thinwall-tubes

have to be filled completely for given max. speed.
For Fixed-Angle-Rotors use caps.

Thickwall-tubes


can be used with or without o-ring alu-caps.
Partially filled without cap,-medium speed.

Thickwall-tubes

with thread have to be filled completely for given
max. speed. Use threaded tube only with screw cap.

Screw Cap

①





Tube & Screw Cap
Assembly

25 46 94/ 4,4ml
25 46 91/ 6,5ml
25 46 88/ 13,5ml
25 46 85/ 38,5ml
25 46 82/ 94ml

O-ring

Alu Cap

②




25 46 89/ 6,5ml thin
25 46 86/ 13,5ml thin
25 46 90/ 13,5ml thin, titan
25 46 87/ 17ml, 13,5ml, thick
25 46 83/ 38,5ml thin
25 46 84/ 38,5ml thick
25 46 92/ 38,5ml thin, titan
25 46 80/ 94ml thin
25 46 81/ 94ml thick


Selfseal

Cap

③



25 47 16/ 2ml
25 47 15/ 4,4ml
25 47 14/ 6,5ml
25 47 13/ 17ml, 13,5ml
25 47 12/ 38,5ml
25 47 11/ 94ml



Tube Removal
Tool

25 49 00

other Tools
on request

Herolab GmbH
Ludwig-Wagner-Str. 12
D-69168 Wiesloch
Tel: +49(0)6222 5802-0
Fax: +49(0)6222 5802-34
E-Mail: herolab_gmbh@t-online.de

HEMOTEC GmbH
Ribistrasse 7
CH-4460 Gelterkinden
Tel: +41(0)61 981 66 33
Fax: +41(0)61 981 40 26
E-Mail: hemotec@bluewin.ch